

**PROJECT:** Medical Center  
**GROUP NO:** 04                      **BATCH:** T5.2306.E1


## A. SYSTEM REQUIREMENTS:

Sr.	Item	Description
1.	Operating system	Windows 7, 8, or higher
2.	Hardware	CPU Cores or better with 2 Gigabytes of RAM or better
3.	Database	SQL Server 2012
4.	Fonts	Standard font Unicode: Arial, Tahoma, etc...
5.	Other requirement	A working internet connection, etc...

## B. Program Installation

### B.1. Setup Ms Sql Server TCP/IP and Port:


1. Open SQL Server Configuration Manager.




2. Choose SQL Server Network Configuration /Protocols for MSSQLSERVER In SQL Server Configuration Manager window.
3. Ensure TCP/IP is enabled.
4. Ensure TCP/IP Port is set to 1433.

---

(\*) - Installation documentation is included with the installation disc, should follow the proper format which have clear and complete content for anyone to perform.  
- The Center performs the installation and give a notification of the protection list for successfully. Project not successfully installed is Fail.




##### 5. Restart MSSQLSERVER in SQL Server Services.




- (\*)
- Installation documentation is included with the installation disc, should follow the proper format which have clear and complete content for anyone to perform.
  - The Center performs the installation and give a notification of the protection list for successfully. Project not successfully installed is Fail.

## B.2. Generated database from script.sql

### 1. Open script.sql




### 2. Create new database




- (\*)
- Installation documentation is included with the installation disc, should follow the proper format which have clear and complete content for anyone to perform.
  - The Center performs the installation and give a notification of the protection list for successfully. Project not successfully installed is Fail.

### 3. Execute code in **script.sql**



### B.3. Setup Application

1. Copy file **MedicalCenter.rar** to HDD.
2. Extract file to C:\.
3. Ensure the folder structure is as the same as the below image:




4. Run file **run.bat**.

### C. Installation successful

- The login window is displayed.
- Login information:
  - user name: **admin**
  - password: **123123**

---

(\*) - Installation documentation is included with the installation disc, should follow the proper format which have clear and complete content for anyone to perform.  
- The Center performs the installation and give a notification of the protection list for successfully. Project not successfully installed is Fail.



- 
- (\*)
- Installation documentation is included with the installation disc, should follow the proper format which have clear and complete content for anyone to perform.
  - The Center performs the installation and give a notification of the protection list for successfully. Project not successfully installed is Fail.