

PROJECT: Medical Center
GROUP NO: 04 **BATCH:** T5.2306.E1


A. SYSTEM REQUIREMENTS:

Sr.	Item	Description
1.	Operating system	Windows 7, 8, or higher
2.	Hardware	CPU Cores or better with 2 Gigabytes of RAM or better
3.	Database	SQL Server 2012
4.	Fonts	Standard font Unicode: Arial, Tahoma, etc...
5.	Other requirement	A working internet connection, etc...

B. Program Installation


B.1. Setup Ms Sql Server TCP/IP and Port:

1. Open SQL Server Configuration Manager.




2. Choose SQL Server Network Configuration /Protocols for MSSQLSERVER In SQL Server Configuration Manager window.
3. Ensure TCP/IP is enabled.
4. Ensure TCP/IP Port is set to 1433.

(*) - Installation documentation is included with the installation disc, should follow the proper format which have clear and complete content for anyone to perform.
- The Center performs the installation and give a notification of the protection list for successfully. Project not successfully installed is Fail.




5. Restart MSSQLSERVER in SQL Server Services.




- (*)
- Installation documentation is included with the installation disc, should follow the proper format which have clear and complete content for anyone to perform.
 - The Center performs the installation and give a notification of the protection list for successfully. Project not successfully installed is Fail.

B.2. Generated database from script.sql

1. Open script.sql



2. Execute code in script.sql




B.3. Setup Application

1. Copy file **MedicalCenter.rar** to **HDD**.
2. Extract file to **C:**.

(*) - Installation documentation is included with the installation disc, should follow the proper format which have clear and complete content for anyone to perform.
- The Center performs the installation and give a notification of the protection list for successfully. Project not successfully installed is Fail.


3. Ensure the folder structure is as the same as the below image:



4. Run file **run.bat**.

C. Installation successful

- The login window is displayed.
- Login information:
 - user name: **admin**
 - password: **123123**



-
- (*)
- Installation documentation is included with the installation disc, should follow the proper format which have clear and complete content for anyone to perform.
 - The Center performs the installation and give a notification of the protection list for successfully. Project not successfully installed is Fail.

(*) - *Installation documentation is included with the installation disc, should follow the proper format which have clear and complete content for anyone to perform.*
- *The Center performs the installation and give a notification of the protection list for successfully. Project not successfully installed is Fail.*